



**Testing &
Certification
Australia**

14 Nelson Street Chatswood NSW 2067 Australia

Telephone 9410 5111 Facsimile 9415 1567
National 02+, International 612+

Mail: P.O. Box 841 Artarmon NSW 2064 Australia
Email: med@energy.com.au

Test Report No: 76945

EQUIPMENT DETAILS

Type: **Blood Pressure Analysis System**
Manufacturer/Trade Mark: **PWV Medical Pty Ltd**
Model: **BPAS-1/MM2**
Country of Manufacturer: **Australia**
Serial Number: **66-0001-1998**
Rated Input: **220-240V or 100-110V 12VA 50/60Hz**
Classifications: **Class II, Type CF**

SUBMITTER DETAILS

**PWV Medical Pty Ltd
31 Hope Street
Ermington NSW 2115
Australia**

Phone: **(612) 9874 8761**
Fax: **(612) 9874 9022**

STANDARDS APPLIED

IEC601-1:1988, A1:1991 & A2:1995
EN60601-1:1990, A1:1993, A12:1993, A2:1995 & A13:1996

RESULT


The sample was found to be in compliance, except at the client's request the requirements of ISO10993-1 (Biocompatibility), as referred to in Clause 48 was not applied.

Assessed By:


Peter Selvey

Date/s: **26 & 27 Nov. 1998**

Reviewed By:


Scott Wedesweiler

Date: **8 January 1999**